

# RNA Miniprep using 2-Butoxyethanol

by A. Untergasser (contact address and download at [www.untergasser.de/lab](http://www.untergasser.de/lab))

Version: 1.0 - [Print Version \(.PDF\)](#)

**ATTENTION:** This is a low priced protocol. Use it preferably!  
It works well, but it takes three days and involves many toxic chemicals.

1. Add 1% beta-mercaptoethanol to the extraction buffer
2. Grind max. 125 mg tissue into a fine powder
3. Transfer grinded powder into a eppi
4. Add **1 ml of extraction buffer** and vortex well  
Careful, pressure may build up!
5. Spin for **10 min** at top speed
6. Transfer top phase to a fresh eppi
7. Add **500 µl chloroform** and vortex well
8. Spin for **5 min** at top speed
9. Transfer top phase to a fresh eppi
10. Add **500 µl phenol / chloroform** and vortex well
11. Spin for **10 min** at top speed
12. Transfer top phase to a fresh eppi
13. Add **500 µl chloroform** and vortex well
14. Spin for **5 min** at top speed
15. Transfer top phase to a fresh eppi
16. Add **40 µl 3 M NaAc and 1 ml 96% EtOH** and vortex well
17. Store at -80 °C until frozen
18. Spin for **5 min** at top speed at 4 °C
19. Remove supernatant and add **500 µl 70% EtOH** to wash
20. Spin for **5 min** at top speed
21. Remove supernatant and add **500 µl RNA resuspension buffer**
22. Incubate at 65 °C to dissolve pellet
23. Add **200 µl 2-Butoxyethanol** and mix well
24. Incubate on ice for at least 3 hours.
25. Spin for **10 min** at max. speed at 4 °C
26. Transfer supernatant to a fresh eppi
27. Add **300 µl 2-Butoxyethanol** and mix well
28. Incubate on ice over night.

29. Spin for **10 min** at max. speed at **4 °C**
30. Discard the supernatant and add **500 µl 70% EtOH**
31. Spin for **5 min** at max. speed at **4 °C**
32. Discard the supernatant and let pellet briefly dry
33. Dissolve in **1 ml water** (and transfer to an eppi)
34. Precipitate RNA by adding **335 µl 8 M LiCl**
35. Incubate at **4 °C** over night
- .
36. Spin for **30 min** at max. speed at **4 °C**
37. Discard the supernatant and add **100 µl 70% EtOH**
38. Spin for **5 min** at max. speed
39. Discard the supernatant and dissolve in **200 µl water**
40. Add **20 µl 3 M NaAcetate (pH 5.2)** and **500 µl 96% ethanol**
41. Spin for **5 min** at max. speed
42. Store at **-80 °C** until frozen
43. Spin for **30 min** at top speed at **4 °C**
44. Remove supernatant and add **100 µl 70% EtOH** to wash
45. Spin for **5 min** at top speed
46. Remove supernatant and redissolve **50 µl water**
47. Incubate at **65 °C** to dissolve pellet

## **Buffers and Solutions**

### **RNA Extraction Buffer**

Prepare in a 1 liter bottle:

Triisopropylnaphtalene Sulfonic Acid	5.56 g	1%
p-Aminosalicylic Acid	30 g	6%
1 M Tris-HCl (pH 8.0)	50 ml	100 mM
0.5 M EGTA (pH 8.0)	50 ml	50 mM
NaCl	2.9 g	100 mM
SDS	5 g	1%
Polyvinylpyrrolidone (PVP) (M.W. 40,000)	5 g	1%
Polyvinylpolypyrrolidone (PVPP)	15 g	3%
add water to 500 ml		

Autoclave for 15 min, then add:

chloroform	500 ml	
------------	--------	--

The extraction buffer is a mess. The PVPP will not dissolve, the chloroform will always separate and it has a brown color. The best is to put it on a magnetic stirrer and use it from there. It has a cacao color and can be pipetted with 1 ml pipettes were the tip was cut off.

Always add beta-mercaptoethanol fresh.

Before use add Mercaptoethanol:

for 50 ml Buffer add 250  $\mu$ l Mercaptoethanol  
for 30 ml Buffer add 150  $\mu$ l Mercaptoethanol  
for 10 ml Buffer add 50  $\mu$ l Mercaptoethanol  
for 5 ml Buffer add 25  $\mu$ l Mercaptoethanol  
for 1 ml Buffer add 5  $\mu$ l Mercaptoethanol

## **RNA Resuspension Buffer**

Boric Acid	0.78 g	25 mM
1 M Tris-HCl (pH 7.6)	25 ml	50 mM
0.5 M EGTA (pH 8.0)	1.3 ml	1.25 mM
NaCl	2.9 g	100 mM
add water to 500 ml		

## **Required Solutions and Chemicals**

- 3 M Natriumacetate (NaAc) solution (pH 5.2)
- 2-Butoxyethanol (2BE)
- 8 M Lithiumchloride (LiCl)

## **Commented Protocol:**

### **1. Add 1% beta-mercaptoethanol to the extraction buffer**

Because the extraction buffer always separates it is the best to put it on a magnetic stirrer and use it from there. It can be pipetted with 1 ml pipettes where the tip was cut off. Always add beta-mercaptoethanol fresh.

### **2. Grind max. 125 mg tissue into a fine powder**

This is best done in a mortar for big amounts with liquid nitrogen or in a cap-shaker for smaller amounts.

### **3. Transfer grinded powder into a eppi**

### **4. Add 1 ml of extraction buffer and vortex well**

Careful, pressure may build up!

### **5. Spin for 10 min at top speed**

The chloroform is already in the extraction buffer.

### **6. Transfer top phase to a fresh eppi**

### **7. Add 500 µl chloroform and vortex well**

### **8. Spin for 5 min at top speed**

**9. Transfer top phase to a fresh eppi**

**10. Add 500 µl phenol / chloroform and vortex well**

Just mix phenol and chloroform 1:1 and mix well before use. It also separates after some time.

**11. Spin for 10 min at top speed**

**12. Transfer top phase to a fresh eppi**

**13. Add 500 µl chloroform and vortex well**

**14. Spin for 5 min at top speed**

**15. Transfer top phase to a fresh eppi**

**16. Add 40 µl 3 M NaAc and 1 ml 96% EtOH and vortex well**

This step precipitates all nucleic acids. The pellet has a brown color.

**17. Store at -80 °C until frozen**

**18. Spin for 5 min at top speed at 4 °C**

**19. Remove supernatant and add 500 µl 70% EtOH to wash**

**20. Spin for 5 min at top speed**

**21. Remove supernatant and add 500 µl RNA resuspension buffer**

**22. Incubate at 65 °C to dissolve pellet**

Pipetting up and down can help as well.

**23. Add 200 µl 2-Butoxyethanol and mix well**

In this first step we precipitate polysaccharides at a low 2-BE concentration.

**24. Incubate on ice for at least 3 hours.**

The publication incubates for more than 30 min, it seems both is working.

**25. Spin for 10 min at max. speed at 4 °C**

It will form a brown pellet containing the sugars.

**26. Transfer supernatant to a fresh eppi**

**27. Add 300 µl 2-Butoxyethanol and mix well**

In this second step we precipitate nucleotides at a high 2-BE concentration.

**28. Incubate on ice over night.**

The publication incubates for more than 30 min and I extended it up to 24 hours, it seems as if both is working.

**29. Spin for 10 min at max. speed at 4 °C**

It will form a brown pellet containing the nucleotides.

**30. Discard the supernatant and add 500 µl 70% EtOH**

**31. Spin for 5 min at max. speed at 4 °C**

The pellet does not look different.

**32. Discard the supernatant and let pellet briefly dry**

Not too long, over dried nucleotides are difficult to dissolve in water.

**33. Dissolve in 1 ml water (and transfer to an eppi)**

**34. Precipitate RNA by adding 335 µl 8 M LiCl**

**35. Incubate at 4 °C over night**

Lithiumchlorid precipitates RNA, the timing seems to be critical as well. The majority should be precipitated after four hours.

**36. Spin for 30 min at max. speed at 4 °C**

It will form a brown pellet containing the nucleotides.

**37. Discard the supernatant and add 100 µl 70% EtOH**

**38. Spin for 5 min at max. speed**

It will form a brown pellet containing the nucleotides.

**39. Discard the supernatant and dissolve in 200 µl water**

**40. Add 20 µl 3 M NaAcetate (pH 5.2) and 500 µl 96% ethanol**

**41. Spin for 5 min at max. speed**

It will form a brown pellet containing the nucleotides.

**42. Store at -80 °C until frozen**

**43. Spin for 30 min at top speed at 4 °C**

**44. Remove supernatant and add 100 µl 70% EtOH to wash**

**45. Spin for 5 min at top speed**

**46. Remove supernatant and redissolve 50 µl water**

**47. Incubate at 65 °C to dissolve pellet**

Pipetting up and down can help as well.

### **Known Issues:**

- This extraction method takes three days and lots of hand on time.
- The buffer is a mess and some compounds are difficult to get.

### **References and Comments:**

This protocol should allow to extract RNA from difficult tissues, which are rich on sugars and phenolic compounds. I got the protocol from J. Schaart and tried it out for strawberry roots. I got great RNA, but still the qPCR had some problems. If you have problems extracting RNA, give it a try, but it is a very labor intensive protocol.

The magic step is the precipitation of sugars with 2-BE. Alternatively you could use your own extraction buffer and just add the 2-BE steps to your protocol. Good luck!

### **How to cite this page in publications:**

This document can be cited like this:

Untergasser A. "RNA Miniprep using 2-Butoxyethanol" *Untergasser's Lab*. Summer 2008. (include here the date when you accessed these page).

<[http://www.untergasser.de/lab/protocols/miniprep\\_rna\\_butoxyethanol\\_long\\_v1\\_0.htm](http://www.untergasser.de/lab/protocols/miniprep_rna_butoxyethanol_long_v1_0.htm)>.

It is heavily based on the following publication. Maybe you prefer to cite the original paper:

Schultz DJ, Craig R, Cox-Foster DL, Mumma RO, Medford JI.

RNA Isolation from Recalcitrant Plant Tissue.

Plant Molecular Biology Reporter, 310-316, 12 (4) 1994.

### **Please Do Not Reprint This Article:**

This article is copyrighted. Please do not reproduce this article in whole or part, in any form, without [obtaining my written permission](#).